

Nasal Mask

SwingFit™.
Performance defined.



Fisher & Paykel
HEALTHCARE



Discover the definition of performance

We know that for a sound sleep, mask performance is crucial. This is why our industry-leading engineers examined mask design from every angle to create a mask that you and your patients can trust. The result: F&P Nova™ Nasal - it's performance defined.



Combining innovative technologies

Packed with features working in synergy with the intuitive SwingFit headgear, F&P Nova Nasal delivers all-round performance. Discover a mask that's easy to fit and provides a stable seal, with a clear field of view and minimal noise and draft.





Product Codes

Code	Description
NVN1SZ	Nova Nasal Mask – Small
NVN1MZ	Nova Nasal Mask – Medium
NVN1LZ	Nova Nasal Mask – Large

NVN1SMLZ Nova Nasal Mask – Fit pack – Small, Medium, Large

Spare Parts

Code	Description
400NVN111	Nova Nasal Spare Cushion – Small
400NVN112	Nova Nasal Spare Cushion – Medium
400NVN113	Nova Nasal Spare Cushion – Large
400NVN121	Nova Nasal Headgear – Small
400NVN122	Nova Nasal Headgear – Medium-Large
400NVN123	Nova Nasal Back Panel – Small*
400NVN124	Nova Nasal Back Panel – Medium-Large *
400NVN141	Nova Nasal Frame Assembly Spare
400NVN151	Nova Nasal Diffuser Spare
400NVN161	Nova Nasal Tube Clip Spare
400NVN162	Nova Nasal Headgear Pad Spare



Tube clip

Headgear pad

Back panel

Headgear

Diffuser

Cushion

Frame assembly
(frame + tube + swivel)



F&P myMask™ App

Instructs your patients on mask use with step-by-step fitting, fine-tuning and cleaning videos. Download the F&P myMask App. Scan the QR code or download directly.



*HEADGEAR INCLUDES BACK PANEL. BACK PANEL CAN ALSO BE PURCHASED SEPARATELY.

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS TO USE. AS PRESCRIBED OR STOWED. F&P myMask is a prescription only. THIS PRODUCT MAY NOW BE THE SUBJECT OF A FOLLOW-UP STATEMENT: "I HEREBY STATE THAT F&P OF FISHER & PAYKEL HEALTHCARE LIMITED. FOR YOUR HEALTHCARE PROVIDER WILL ADVISE YOU ON THE BEST MASK TO WORK." INTERNAL EVALUATION OF INFORMATION, PLEASE SEE WWW.FPHCARE.COM/IFP. PRODUCT IS SUITABLE FOR YOU/YOUR CONDITIONING. IF YOU HAVE ANY DOUBTS BY FISHER & PAYKEL HEALTHCARE LIMITED. A SERVICE MARK OF APPLE INC., REGISTERED PATENT, TALK TO YOUR HEALTHCARE PROVIDER. 43 PARTICIPANTS IN THE UNITED STATES, 2022. IN OTHER COUNTRIES AND IN THE UNITED STATES, 2022. AND THE GOOGLE PLAY LOGO ARE TRADEMARKS OF GOOGLE LLC. 63666 REV A © 2025 FISHER & PAYKEL HEALTHCARE LIMITED.

Fisher & Paykel
HEALTHCARE