

Mask/Device Compatibility List

English

This document specifies the masks compatible with ResMed's positive airway pressure devices and the mask setting, if one is required. For full details regarding the correct use of these devices, please refer to the specific User Guides.

Note: Not all products are available in all regions.

Full Face Masks 1

Device	AcuCare F1-4	AirFit™ F40 ⁴	AirTouch™ F30i ^{1,4}	AirFit™ F30i ¹ / AirFit™ F30i Non Magnetic ¹	AirFit™ F30	AirTouch™ F20 / AirTouch™ F20 for Her	AirFit™ F20 / AirFit™ F20 for Her/ AirFit™ F20 Non Magnetic	AirFit™ F10 / AirFit™ F10 for Her
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	X	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 11 ASV PaceWave	X	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirSense™ 11	N/V	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirMini™	X	X	X	X	✓	✓	✓	X
AirStart™ 10 / AirSense™ 10	N/V	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	N/V	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 CS Pacewave	N/V	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 CS-A Pacewave	N/V	N/V	N/V	N/V	Full Face	Full Face	Full Face	Full Face
Lumis™	N/V	Pillows ³	Full Face ³	Full Face ³	Full Face ³	Full Face	Full Face	Full Face
S9 AutoSet™ / S9 Elite™	N/V	Pillows	Full Face	Full Face	N/V	Full Face	Full Face	Full Face
S9 Escape™	N/V	Pillows	Full Face	Full Face	N/V	Full Face	Full Face	Full Face
AutoSet CS™ 2	N/V	N/V	N/V	N/V	N/V	N/V	N/V	FULL FACE
AutoSet CS™ (S9)	N/V	Pillows	N/V	N/V	N/V	Full Face	Full Face	Full Face
AutoSet CS™ (S9 PaceWave)	N/V	Pillows	N/V	N/V	N/V	Full Face	Full Face	Full Face
AutoSet CS-A™ (S9 PaceWave)	N/V	N/V	N/V	N/V	N/V	Full Face	Full Face	Full Face
C-Series Tango™	N/V	N/V	N/V	N/V	N/V	N/V	N/V	No setting required. Consult your clinician before changing to a different mask type.
S7™ Lightweight	N/V	N/V	N/V	N/V	N/V	N/V	N/V	No setting required. Consult your clinician before changing to a different mask type.
S7 Elite™ / S7 AutoSet Spirit	N/V	N/V	N/V	N/V	N/V	N/V	N/V	MIR FULL
S7 AutoSet Respond™	N/V	N/V	N/V	N/V	N/V	N/V	N/V	MIRAGE
S8™ Auto 25	N/V	N/V	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™	N/V	N/V	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Vantage™	N/V	N/V	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL
S8 Escape™	N/V	N/V	N/V	N/V	N/V	N/V	N/V	Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.

Device	AcuCare F1-4	AirFit™ F404	AirTouch™ F30i	AirFit™ F30i / AirFit™ F30i Non Magnetic1	AirFit™ F30	AirTouch™ F20 / AirFit™ F20 for Her	AirFit™ F20 / AirFit™ F20 Non Magnetic	AirFit™ F10 / AirFit™ F10 for Her
S8 Escape™ II S8	N/	N/V	N/V	N/V	N/	Requires pressure calibration at 12cm H2O. Consult your clinician and check the user guide.		
Elite™ / S8 Elite™ II S8	V	N/V	N/V	N/V	V	MIR FULL	MIR FULL	MIR FULL
AutoSet Spirit™ II S8	N/	N/V	N/V	N/V	N/	MIR FULL	MIR FULL	MIR FULL
Compact™ VPAP™ III / VPAP™ III ST VPAP™	V	N/V	N/V	N/V	V	N/V	N/V	No setting required. Consult your clinician before changing to a different mask
III ST-A VPAP™ IV / VPAP™ IV ST VPAP™	N/	N/V	N/V	N/V	N/	N/V	N/V	
III ST-A VPAP™ IV / VPAP™ IV ST VPAP™	V	N/V	N/V	N/V	V	N/V	N/V	MIR FULL
VPAP™ IV ST VPAP™	N/	N/V	N/V	N/V	N/	MIR FULL	MIR FULL	MIR FULL
Adapt SV VPAP™ Adapt	V	N/V	N/V	N/V	V	N/V	N/V	FULL FACE
(S9) VPAP™ Adapt (S9	N/	Pillows	Full Face	Full Face	N/	Full Face	Full Face	Full Face
PaceWave) VPAP™	V	Pillows	N/V	Full Face	V	Full Face	Full Face	Full Face
Malibu VPAP™ Auto 25	N/	N/V	N/V	N/V	N/	N/V	N/V	FULL FACE
Auto 25 (S9) VPAP™ S / VPAP™ ST VPAP™ S	V	N/V	N/V	N/V	V	MIR FULL	MIR FULL	MIR FULL
VPAP™ ST VPAP™ S	N/	Pillows	N/V	Full Face	N/	Full Face	Full Face	Full Face
(S9) / VPAP™ ST (S9)	V	N/V	N/V	N/V	V	MIR FULL	MIR FULL	MIR FULL
VPAP™ ST-A (S9)	N/	Pillows	Full Face	Full Face	N/	Full Face	Full Face	Full Face
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	V	N/V	N/V	N/V	V	Full Face	Full Face	Full Face
VPAP™ ST-A (S9 iVAPS)	N/	N/V	N/V	N/V	N/	Full Face	Full Face	Full Face
S9 VPAP™ Tx	V				V			
VPAP™ Tx	N/	Pillows2	Full Face2	Full Face2	N/	Full Face	Full Face	Full Face
Stellar™	V	N/V	N/V	N/V	V	N/V	N/V	FULL FACE
VS Ultra™ / VS	N/	N/V	N/V	N/V	N/	FULL FACE	FULL FACE	FULL FACE
Integra™	V	N/V	N/V	N/V	V	N/V	N/V	
VS III™	N/	N/V	N/V	N/V	N/	N/V	N/V	✓
Astral™	N/V	N/V	N/V	N/V	V	FULL FACE	FULL FACE	✓
	FULL FACE				N/			N/V

¹ SmartStop may not operate as described when using this mask with some CPAP or bi-level devices.

² Excluding PAC and iVAPS modes. Only compatible for pressure range 4 - 25 cm H2O.

³ This combination of mask and device has been determined safe for use based on the evaluation of accuracy of monitoring of pressure, volume, minute ventilation and respiratory rate, and its alarms if applicable. Small variations in performance might be observed, As per expected clinical practice, the effectiveness of ventilation and alarms should be verified by the clinician.

⁴ SmartStart may not operate as described when using this mask with some CPAP or bi-level devices.

N/V = The mask has not been validated for this device.

✓ = Compatible X = Not Compatible

Full Face Masks 2

Device	N/	Mirage™	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ FX / Quattro™ FX for Her	Mirage Liberty™	Hospital Full Face
	V			V				
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 STV	N/	N/V	N/V	Full Face	N/	Full Face	Full Face	Full Face
AirCurve™ 11 ASV PaceWave	N/	N/V	N/V	Full Face	N/	Full Face	Full Face	Full Face
AirSense™ 11	V	N/V	N/V	Full Face	V	Full Face	Full Face	Full Face
AirMini™		X	X	X	N/	X	X	X

Device	Mirage™	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ Quattro™	FX / Quattro™ Air / Quattro™	FX for Her Air for Her	Mirage Liberty™	Hospital Full Face
AirStart™ 10 / AirSense™ 10	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
AirCurve™ 10 CS Pacewave	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
AirCurve™ 10 CS-A Pacewave	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
Lumis™	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
S9 AutoSet™ / S9 Elite™	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
S9 Escape™	FULL FACE	FULL FACE	FULL FACE	FULL FACE		FULL FACE		FULL FACE1	FULL FACE
AutoSet CS™ 2	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
AutoSet CS™ (S9)	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
AutoSet CS™ (S9 PaceWave)	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
AutoSet CS-A™ (S9 PaceWave)			No setting required. Consult your clinician before changing to a different mask type.						
C-Series Tango™	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR	MIR
S7™ Lightweight	STANDARD	ULTRA	MIRAGE	MIRAGE		MIRAGE		FULL	FULL
S7 Elite™ / S7 AutoSet Spirit	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR	MIR
S7 AutoSet Respond™	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		FULL	FULL
S8™ Auto 25	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR	MIR
S8 AutoSet Spirit™			Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.					FULL	FULL
S8 AutoSet Vantage™			Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.					MIR	MIR
S8 Escape™	MIR	MIR	MIR	MIR FULL		MIR FULL		FULL	FULL
S8 Escape™ II	FULL	FULL	FULL	MIR FULL		MIR FULL		FULL	FULL
S8 Elite™ / S8 Elite™ II	MIR	MIR	MIR	No setting required. Consult your clinician before changing to a different mask type.					FULL
S8 AutoSet Spirit™ II	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR FULL	MIR FULL
S8 Compact™	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR FULL	MIR FULL
VPAP™ III / VPAP™ III ST	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR FULL	MIR FULL
VPAP™ III ST-A	FULL FACE	FULL FACE	FULL FACE	FULL FACE		FULL FACE		FULL FACE1	FULL FACE
VPAP™ IV / VPAP™ IV ST	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
VPAP™ Adapt SV	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
VPAP™ Adapt (S9)	FULL FACE	FULL FACE	FULL FACE	FULL FACE		FULL FACE		FULL FACE	FULL FACE
VPAP™ Adapt (S9 PaceWave)	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR FULL	MIR FULL
VPAP™ Malibu	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
VPAP™ Auto 25	MIR FULL	MIR FULL	MIR FULL	MIR FULL		MIR FULL		MIR FULL	MIR FULL
Auto 25 (S9)	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
VPAP™ S / VPAP™ ST	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
VPAP™ S (S9) / VPAP™ ST (S9)	N/V	N/V	Full Face	Full Face		Full Face		Full Face	Full Face
VPAP™ ST-A (S9)	N/V	N/V	Full Face					Full Face	Full Face
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS) S9	N/V	N/V	FULL FACE	Full Face		Full Face		FULL FACE	FULL FACE
VPAP™ Tx VPAP™ Tx Stellar™			FULL FACE	FULL FACE		FULL FACE		FULL FACE	FULL FACE

Device	Mirage™	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ Quattro™	FX / Quattro™ Air / Quattro™	FX for Her Air for Her	Mirage Liberty™	Hospital Full Face
VS Ultra™/ VS Integra™	N/	N/	✓	✓		N/V		✓	✓
VS II™	V	V	✓	✓		N/V		✓	✓
Astral™	N/	N/	FULL FACE	FULL FACE		FULL FACE		FULL FACE	FULL FACE
	V	V							
	N/	N/							
	V	V							

Nasal Masks 1

Device	AirTouch™ N30i1	AirFit™ N30i1 for AirMini	AirFit™ N30i1 AirFit™ N30	AirFit™ N20 Classic	AirFit™ N20 / AirFit™ N20 for Her	AirTouch™ N20 / AirFit™ N20 for Her	AirFit™ N10/Swift™ FX Nand	FX Nand	Mirage™ FX / Mirage™ FX™ for Her	Mirage™ Active™ LT	Mirage™ SoftGel	Hospital Nasal
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™	Pillows	X	Pillows	Pillo	Pillo	Pillow	Pillows	Pillows	Nasa	Nasa	X	Nasal
11 ST	Pillows	X	Pillows	ws	ws	s	Pillows	Pillows	I	I	X	Nasal
AirCurve™ 11 ASV PaceWave	Pillows	X	Pillows	Pillo	Pillo	Pillow	Pillows	Pillows	Nasa	Nasa	N/V	Nasal
AirSense™ 11	X	✓	X	✓	✓	✓	X	X	X	X	X	X
AirMini™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirStart™ 10 / AirSense™ 10	Pillows	X	Pillows	Pillo	Pillo	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 VAuto /	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 S /	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 ST	N/V	X	N/V	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 CS Pacewave	Pillows ³	X	Pillows ³	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 CS-A Pacewave	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
Lumis™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 AutoSet™ / S9 Elite™	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V	VISTA1	ACTIVA	ACTIVA	VISTA
S9 Escape™	N/V	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ 2	N/V	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ (S9)	N/V	X	N/V	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ (S9 PaceWave)	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V				
AutoSet CS-A™ (S9	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V	No setting required. Consult your clinician before changing to a			
PaceWave)	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V	different mask type. No setting required. Consult your clinician			
C-Series Tango™	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V	before STANDARD to a STANDARD mask type.	STANDARD	STANDARD	STANDARD
S7™ Lightweight	N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S7 Elite™ / S7 AutoSet Spirit	N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S7 AutoSet Respond™	N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8™ Auto 25	N/V	X	N/V	N/V	N/V							
S8 AutoSet Spirit™	N/V	X	N/V	N/V		N/V	Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.					
S8 AutoSet Vantage™	N/V	X	N/V	N/V			Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.					
S8 Escape™	N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 Escape™ II	N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 Elite™ / S8 Elite™ II	N/V	X	N/V	N/V	N/V	N/V	No setting required. Consult your clinician before changing to a different mask type.					
S8 AutoSet Spirit™ II	N/V	X	N/V	N/V	N/V	N/V	MIRAGE	MIRAGE	STANDARD	STANDARD	STANDARD	N/V
S8 Compact™	N/V	X	N/V	N/V	N/V	N/V	MIRAGE	MIRAGE	STANDARD	STANDARD	STANDARD	N/V
VPAP™ III / VPAP™ III ST		X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
VPAP™ III ST-A					N/V	N/V	N/V	N/V	VISTA1	ACTIVA	ACTIVA	VISTA
VPAP™ IV / VPAP™ IV ST												
VPAP™ Adapt SV												

Device			AirTouch™ N30i ¹	AirFit™ N30 for AirMini	AirFit™ N30i ¹ AirFit™ N30 Classic	AirFit™ N20	AirFit™ N20 / AirFit™ N20 for Her	AirTouch™ N20 / AirFit™ N20 for Her	AirFit™ N10/Swift™ FX Nand	Mirage™ FX / Mirage™ FX™ for Her	Mirage™ LT	Mirage™ SoftGel	Hospital Nasal	
VPAP™	Adapt	(S9)	N/V	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™	Adapt	(S9	N/V	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
PaceWave) VPAP™			N/V	X	N/V	N/V	N/V	N/V	SWIFT	SWIFT	VISTA1	N/V	N/V	MERIDIAN
Malibu VPAP™ Auto 25			N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
Auto 25 (S9) VPAP™ S /			N/V	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST VPAP™ S			N/V	X	N/V	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
(S9) / VPAP™ ST (S9)			Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST-A (S9)			N/V	X	N/V	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)			N/V	X	N/V	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 VPAP™ Tx			Pillows2	X	Pillows2	N/V	Pillows	N/V	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Tx			N/V	X	N/V	N/V	N/V	N/V	PILLOWS	PILLOWS	NASAL	NASAL	NASAL	NASAL
Stellar™			N/V	X	N/V	PILLOWS	PILLOWS	PILLOWS	PILLOWS1	PILLOWS1	NASAL	NASAL	NASAL	NASAL
VS Ultra™ / VS			N/V	X	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V
Integra™			N/V	X	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V
VS III™			N/V	X	N/V	Pillows	Pillows	Pillows	N/V	N/V	Nasal	✓	Nasal	✓
Astral™												N/V		Nasal

¹ SmartStop may not operate as described when using this mask with some CPAP or bi-level devices.

² Excluding PAC and iVAPS modes. Only compatible for pressure range 4 - 25 cm H2O.

³ This combination of mask and device has been determined safe for use based on the evaluation of accuracy of monitoring of pressure, volume, minute ventilation and respiratory rate, and its alarms if applicable. Small variations in performance might be observed, As per expected clinical practice, the effectiveness of ventilation and alarms should be verified by the clinician.

N/V = The mask has not been validated for this device.

✓ = CompatibleX = Not Compatible

Nasal Masks 2

Device	Pixi™	Mirage Kidsta™	Mirage Micro™	Mirage Vista™	Ultra Mirage II	Mirage	Mirage Activa™	Modular	Papillon™	Meridian
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	Pediatric	X	X	X	X	N/V	X	N/V	N/V	X
AirCurve™ 11 ASV PaceWave		X	X	X	X	N/V	X	N/V	N/V	X
AirSense™ 11	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V
AirMini™	N/V	X	X	X	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	X	Nasal	Nasal	Nasal	Nasal	N/V	Nasal	N/V	N/V	Nasal
AirCurve™ 10 VAuto / AirCurve™ 10 S /	N/V	Nasal	Nasal	Nasal	Nasal	N/V	Nasal	N/V	N/V	Nasal
AirCurve™ 10 ST	N/V									
AirCurve™ 10 CS Pacewave		Nasal	Nasal	Nasal	Nasal	N/V	Nasal	N/V	N/V	Nasal
AirCurve™ 10 CS-A Pacewave	N/V	Nasal	Nasal	Nasal	Nasal	N/V	Nasal	N/V	N/V	Nasal
Lumis™	N/V	Nasal	Nasal	Nasal	Nasal	N/V	Nasal	N/V	N/V	Nasal
S9 AutoSet™ / S9 Elite™	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
S9 Escape™	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ 2	N/V	N/V	VISTA1	VISTA	ULTRA	N/V	ACTIVA	N/V	PAPILLON	VISTA
AutoSet CS™ (S9)	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal

Device	Pixi™	Mirage Kidsta™	Mirage Micro™	Mirage Vista™	Ultra Mirage II	Mirage	Mirage Activa™	Modular	Papillon™	Meridian
AutoSet CS™ (S9 PaceWave) / AutoSet CS-A™ (S9 PaceWave)	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
C-Series Tango™	N/V	N/V	N/V		No setting required. Consult your clinician before changing to a different mask type.					N/V
S7™ Lightweight	N/V	N/V	STANDARD		No setting required. Consult your clinician before changing to a different mask type.					N/V
S7 Elite™ / S7 AutoSet	N/V	N/V	STANDARD	STANDA	ULTRA	MIRAGE	STANDARD	STANDARD	STANDARD	N/V
Spirit	N/V	N/V	STANDARD	RD	ULTRA	MIRAGE	STANDARD	STANDARD	STANDARD	STANDARD
S7 AutoSet Respond™	N/V	N/V	STANDARD	STANDA	ULTRA	MIRAGE	ACTIVA	N/V	N/V	STANDARD
S8™ Auto 25	N/V	N/V	STANDARD	RD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 AutoSet Spirit™	N/V	N/V		STANDA	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	
S8 AutoSet Vantage™	N/V	N/V		RD	Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.					
S8 Escape™	N/V	N/V		STANDA	Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.					
S8 Escape™ II	N/V	N/V	STANDA	RD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDA	STANDA
S8 Elite™ / S8 Elite™ II	N/V	N/V	RD	RD	ULTRA	MIRAGE	ACTIVA	STANDARD	RD	RD
S8 AutoSet Spirit™ II	N/V	N/V	STANDA	RD	No setting required. Consult your clinician before changing to a different mask type.				STANDA	STANDA
S8 Compact™	N/V	STANDARD	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	RD/N/V	RD/N/V
VPAP™ III / VPAP™ III ST	N/V	N/V	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	N/V	N/V
VPAP™ III ST-A	N/V	N/V	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/V	N/V	STANDARD
VPAP™ IV / VPAP™ IV ST	N/V	N/V	VISTA1	VISTA	ULTRA	N/V	ACTIVA	N/V	N/V	VISTA
VPAP™ Adapt SV	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
VPAP™ Adapt (S9)	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
VPAP™ Adapt (S9 PaceWave)	N/V	N/V	MERIDIAN	VISTA	ULTRA	N/V	N/V	N/V	N/V	MERIDIAN
VPAP™ Malibu	N/V	N/V	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/V	N/V	STANDARD
VPAP™ Auto 25	N/V	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	STANDARD
Auto 25 (S9)	N/V	N/V	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/V	N/V	STANDARD
VPAP™ S / VPAP™ ST	Pediatric		Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
VPAP™ S (S9) / VPAP™										
ST (S9) / VPAP™ ST-A (S9)	Pediatric	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	Nasal
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	Pediatric	N/V	Nasal	Nasal	Nasal Ultra	N/V	Nasal	N/V	N/V	Nasal
S9 VPAP™ Tx	Pediatric	NASAL	NASAL	NASAL	NASAL-ULTRA	N/V	NASAL	N/V	N/V	NASAL
VPAP™ Tx	N/V	N/V	N/V	N/V	ULTRA	N/V	NASAL	N/V	N/V	N/V
Stellar™	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V
VS Ultra™ / VS Integra™	Pediatric	Nasal	Nasal	Nasal	Nasal Ultra	N/V	N/V	N/V	N/V	N/V

VS III™
 N/V = The mask has not been validated for this device.
 Astral™
 ✓ = Compatible X = Not Compatible

Nasal Pillows

Device	AirFit™ P30i1	AirFit™ P10 For AirMini™	AirFit™ P10 / AirFit™ P10 for Her	Mirage Swift™	Mirage Swift™	Swift™ FX / Swift™ FX for Her / Swift™ FX Bella / Swift™ FX Bella Gray	Swift™ LT	Swift™ LT for Her/ Swift™ LT-F
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	Pillows	X	Pillows2	N/V	N/V	Pillows	Pillows	Pillows
AirCurve™ 11 ASV PaceWave	Pillows	X	Pillows	N/V	N/V	Pillows	Pillows	Pillows
AirSense™ 11	X	✓	X	X	X	X	X	X
AirMini™	Pillows	X	Pillows2	N/V	Pillows	Pillows	Pillows	Pillows
AirStart™ 10 / AirSense™ 10	Pillows	X	Pillows2	N/V	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 VAuto /	Pillows	X	Pillows2	N/V	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 S /	Pillows	X	Pillows2	N/V	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 ST	N/V	X	Pillows2	N/V	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 CS Pacewave	Pillows4	X	Pillows	N/V	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 CS-A Pacewave	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows
Lumis™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows
S9 AutoSet™ / S9 Elite™	N/V	X	N/V	Pillows	Pillows	Pillows	Pillows	N/V
S9 Escape™	Pillows	X	Pillows2	N/V	N/V	N/V	N/V	Pillows
AutoSet CS™ 2	Pillows	X	Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows
AutoSet CS™ (S9)	N/V	X	Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows
AutoSet CS™ (S9 PaceWave)	N/V	X		Pillows	Pillows	Pillows	Pillows	
AutoSet CS-A™ (S9)	N/V	X		No setting required. Consult your clinician before changing to a different mask type.				
PaceWave)	N/V	X		No setting required. Consult your clinician before changing to a different mask type.				
C-Series Tango™	N/V	X	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
S7™ Lightweight	N/V	X	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
S7 Elite™ / S7 AutoSet Spirit	N/V	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S7 AutoSet Respond™	N/V	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8™ Auto 25	N/V	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8 AutoSet Spirit™	N/V	X		Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.				
S8 AutoSet Vantage™	N/V	X		Requires pressure calibration at 12 cm H2O. Consult your clinician and check the user guide.				
S8 Escape™	N/V	X	SWI	SWI	SWIFT	SWIFT	SWIFT	SWI
S8 Escape™ II	N/V	X	FT	FT	SWIFT	SWIFT	SWIFT	FT
S8 Elite™ / S8 Elite™ II	N/V	X	SWI	SWI	No setting required. Consult your clinician before changing to a different mask type.			SWI
S8 AutoSet Spirit™ II	N/V	X	MIRAGE 2	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
S8 Compact™	N/V	X	MIRAGE 2	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
VPAP™ III / VPAP™ III ST	Pillows	X	SWIFT2	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ III ST-A	Pillows	X	N/V	N/V	N/V	N/V	N/V	N/V
VPAP™ IV / VPAP™ IV ST	N/V	X	Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Adapt SV	N/V	X	Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Adapt (S9)	Pillows	X	SWIFT2	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ Adapt (S9 PaceWave)	N/V	X	SWIFT2	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ Malibu	Pillows	X	Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Auto 25		X	SWIFT2	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
Auto 25 (S9)			Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows

VPAP™ S / VPAP™ ST

VPAP™ S (S9) / VPAP™ ST (S9)

Device	AirFit™ P30i1	AirFit™ P10 For AirMini™	AirFit™ P10 / AirFit™ P10 for Her	Mirage Swift™	Mirage Swift™ LT	Swift™ FX / Swift™ FX for Her / Swift™ FX Bella / Swift™ FX Bella Gray	Swift™ LT	Swift™ LT for Her/ Swift™ LT-F
VPAP™ ST-A (S9)	N/V	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	N/V	X	Pillows2	Pillows	Pillows	Pillows	Pillows	Pillows
S9 VPAP™ Tx	N/V	X	PILLOWS2	PILLOWS	PILLOWS	PILLOWS	PILLOWS	PILLOWS
VPAP™ Tx	N/V	X	PILLOWS1, 2	PILLOWS1	PILLOWS1	PILLOWS1	PILLOWS1	PILLOWS1
Stellar™	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V
VS Ultra™ / VS Integra™	N/V	X	N/V	N/V	N/V	N/V	N/V	N/V
VS III™	N/V	X	N/V	Pillows	Pillows	Pillows	Pillows	Pillows

¹SmartStop may not operate as described when using this mask with some CPAP or bi-level devices.

²Accuracy of tidal volume/minute ventilation may reduce using the AirFit P10 mask in bilevel modes, particularly at higher pressure supports.

³Excluding PAC and iVAPS modes. Only compatible for pressure range 4 - 25 cm H2O.

⁴This combination of mask and device has been determined safe for use based on the evaluation of accuracy of monitoring of pressure, volume, minute ventilation and respiratory rate, and its alarms if applicable. Small variations in performance might be observed, As per expected clinical practice, the effectiveness of ventilation and alarms should be verified by the clinician.

N/V = The mask has not been validated for this device.

✓ = CompatibleX = Not Compatible